



49 Changi South Avenue 2
Singapore 486056
www.medtronic.com

tel +65 6436 5000
fax +65 6776 6355

URGENT: Medical Device Voluntary Removal
Covidien Parietex™ Composite Parastomal Mesh

15 October 2018

Attention: Risk Management Director and O.R. Materials Management
CC: The Chairman Medical Board and relevant Head of Departments

Dear Valued Customer:

The purpose of this letter is to advise you that Medtronic is issuing a voluntary removal for two item codes of its

Covidien Parietex™ Composite Parastomal Mesh.

Medtronic is issuing this voluntary removal following receipt of reports of Parietex™ composite parastomal mesh failure identified several years following parastomal hernia repair using the modified Sugarbaker repair technique. In these reports, Parietex™ composite parastomal mesh failure led to hernia recurrence requiring additional surgical treatment. Symptoms of hernia recurrence may include discomfort, localized pain-free or painful bulging, and possible changes in the overlying skin. Medtronic has received, worldwide, a total of ten reports of mesh failure following use of Parietex™ composite parastomal mesh in the last five years. Patients who have received a Parietex™ composite parastomal mesh for the treatment of a parastomal hernia should continue to receive ongoing monitoring by their healthcare providers for the recurrence of a parastomal hernia.

This voluntary removal is in relation to all lots of the item code listed below. No other item codes of Parietex™ mesh are affected by this action.

| Item Code | Description | Affected lots |
|-----------|---|---------------------------------|
| PCOPM15 | Parietex™ Composite Parastomal Mesh 15 cm | See lots listed on Attachment A |
| PCOPM20 | Parietex™ Composite Parastomal Mesh 20 cm | |

Medtronic requests that you quarantine and return any unused products of the item codes detailed above. Unused products from the affected item codes should be returned as described in the Required Actions section below. If you have distributed Covidien Parietex™ composite parastomal mesh listed above, please promptly forward the information from this letter to those recipients. All unused products from the affected item codes must be returned.

Required Actions:

1. Please quarantine and discontinue use of the affected item codes and lots listed.
2. Please complete the attached Recalled Product Return Form in its entirety. If you do not have any units from the affected lots in your inventory, simply tick “No affected inventory to be returned” and revert back to your local Medtronic Representative.
3. Please return affected product as follows:

- **CUSTOMERS WHO PURCHASED PRODUCT DIRECTLY FROM MEDTRONIC**

Please complete the **Covidien Parietex™ Composite Parastomal Mesh** Recalled Return Form (attached) and return to your local Medtronic representative together with affected units if any.

- **CUSTOMERS WHO PURCHASED PRODUCT FROM A DISTRIBUTOR**


Please complete the **Covidien Parietex™ Composite Parastomal Mesh** Recalled Return Form (attached) and contact your Distributor directly. All affected product and Recalled Product Return Form must be returned through the Distributor.

This notification is being issued or will be notified to relevant regulatory bodies according to applicable regulations. Please communicate this important information within your facility and or other facilities as required. We request that you contact Medtronic if you experienced quality problems or adverse events.

If you have any questions or concerns regarding this recall, please do not hesitate to contact your local Medtronic representative.

We appreciate your attention to this matter and apologize for any inconvenience this issue may have caused.

Sincerely,



Diana Teo
Quality Management System Manager
Medtronic

Attachment A

| Item Code | Description | Affected Lots | | | | |
|-----------|--|---------------|--|----------|----------|----------|
| PCOPM15 | Parietex™ Composite Parastomal Mesh 15cm | PNI0064 | POJ0880X | PPL0082X | PRA2444X | PRI1409X |
| | | PNJ0989 | POK0364X | PQA0483X | PRB1650X | PRJ0480X |
| | | PNK0641 | POK0789X | PQC0100X | PRB2015X | PRK0758X |
| | | PNL0039 | POL0082X | PQD0395X | PRC0392X | PRK1140X |
| | | POA0469X | PPA0226X | PQE0802X | PRC1062X | PRL0129X |
| | | POB0030X | PPA0508X | PQG1150X | PRD0248X | PRL0528X |
| | | POC0036X | PPB0779X | PQH0608X | PRD0547X | PSA0833X |
| | | POD0036X | PPD0409X | PQJ1246X | PRD1170X | PSA1209X |
| | | POE0262X | PPF0181X | PQJ0971X | PRE0479X | PSB0936X |
| | | POF0295X | PPG0035X | PQK0314X | PRF0408X | PSC0186X |
| | | POG0047X | PPG0723X | PQL0150X | PRG1020X | PSD0420X |
| | | PRG0380X | PPI0630X | PQL0479X | PRH0323X | PSE0906X |
| | | POH0069X | PPI1144X | PRA1193X | PRI0530X | PSF0217X |
| | | POI0260X | PPJ0234X | PRA1655X | PRI1166X | PSG0778X |
| | | POI0489X | PPK0523X | | | |
| | | PCOPM20 | Parietex™ Composite Parastomal Mesh 20cm | PNI0065 | POJ0418X | PPJ0698X |
| PNI0402 | POK0366X | | | PPK0524X | PRB2017X | PRK0760X |
| PNJ0496 | POK0790X | | | PPL0353X | PRC0394X | PRK1142X |
| PNJ0990 | PPA0227X | | | PQC0102X | PRC1064X | PRL0131X |
| PNL0037 | PPA0509X | | | PQD0397X | PRD0250X | PRL0530X |
| PNL0631 | PPB0780X | | | PQE0187X | PRD0549X | PSA0835X |
| POA0963X | PPC0578X | | | PQG0391X | PRD1172X | PSA1207X |
| POB0797X | PPD0410X | | | PQG1152X | PRE0481X | PSB0938X |
| POC0726X | PPE0215X | | | PQH0610X | PRE1319X | PSC0188X |
| POD0033X | PPE0509X | | | PQJ1248X | PRF0410X | PSD0422X |
| POD1016X | PPF0179X | | | PQJ0973X | PRG0382X | PSE0908X |
| POE0149X | PPF0836X | | | PQK0316X | PRG1022X | PSF0219X |
| POE0727X | PPH0355X | | | PQL0481X | PRH0325X | PSG0780X |
| POF0298X | PPI0631X | | | PRA1195X | PRI0531X | RQA1121X |
| POH0070X | PPI1145X | | | PRA1657X | PRI1168X | |
| POI0261X | PPJ0235X | | | PRA2446X | PRI1411X | |

Attachment B

COVIDIEN™
Parietex™
Composite Parastomal Mesh
Monofilament Polyester with Absorbable Collagen Film

15 cm PCOPM15

CE

| | |
|---|---|
| COVIDIEN™ Parietex™ XXXXXXXXXX Use by YYYY-MM-DD | COVIDIEN™ Parietex™ XXXXXXXXXX Use by YYYY-MM-DD |
| COVIDIEN™ Parietex™ XXXXXXXXXX Use by YYYY-MM-DD | COVIDIEN™ Parietex™ XXXXXXXXXX Use by YYYY-MM-DD |

XXXXXXXXXX Use by YYYY-MM-DD

XXXXXXXXXX
XXXXXXXXXX
XXXXXXXXXX

Item code

Lot number

COVIDIEN™ Parietex™ Composite Parastomal Mesh Monofilament Polyester with Absorbable Collagen Film



49 Changi South Avenue 2
Singapore 486056
www.medtronic.com

tel +65 6436 5000
fax +65 6776 6355

RECALLED PRODUCT RETURN FORM
Covidien Parietex™ Composite Parastomal Mesh

ALL CUSTOMERS PLEASE COMPLETE THE FORM IN ITS ENTIRETY

| Customer Contact Details | Medtronic Contact Details |
|--------------------------|---------------------------|
| Hospital / HCP: | By E-mail: |
| Address: | By Post: |
| Telephone no: | Telephone: |
| Fax no: | |
| E-mail: | |

Even if you have no affected stock to returned, please complete this form and table; return to your local Medtronic contacts stated above

Please tick (either one) only:

- No, no affected inventory to be returned.
- Yes, affected inventory to be returned (Please fill up the below table if chosen)

| Item code | Product Description | Lot number | QTY (Each) |
|-----------|---------------------|------------|------------|
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |

I acknowledge receipt of the Covidien Parietex™ composite parastomal mesh removal notification dated 15 October 2018, and understand the instructions provided.

Name: _____ (print) Signature: _____ Stamp: _____ Date: _____