



Biocity Nottingham, Pennyfoot Street, Nottingham NG1 1G, UK  
EU Authorised Representative for Milestone Scientific Inc

### **Field Safety Notice**

Milestone "Wand STA" Hand Piece Part STA-5050-305 Lot # 140103  
2015-05-08

**Type of action:** Examine part and or return for replacement

**Date:** 8 May 2015

**Attention:** Customers of Milestone "Wand STA" Hand Piece Part STA-5050-305 Lot # 140103

**Affected devices:**

The affected devices are Milestone "Wand STA" Hand Piece Part STA-5050-305 Lot # 140103

**Description of the problem:**

A very small proportion of product supplied as Milestone "Wand STA" Hand Piece Part STA-5050-305 Lot # 140103 are in fact Milestone "Wand" handpieces, model WA-2050-305. This is as a result of an incomplete line clearance during a packaging run.

The "Wand" handpiece model WA-2050-305 does not fit the Wand STA system and cannot be used with it, so the problem has no associated patient risk.

**Action to be taken by the user:**

Distributor shall immediately cease distribution of the affected lot and quarantine stock of this lot for checking or replacement.

Distributor can opt to return the affected lot for replacement, or check the stock they hold using a protocol to be supplied by the manufacturer.

Purchasers holding stock from Distributors can return stock for replacement.

Please maintain awareness on this notice and resulting action until all stock has been checked or replaced.




Michael Gamlen  
[michael@pharmdservices.com](mailto:michael@pharmdservices.com) +44 7802 895193

EU Authorised Representative of Milestone Scientific  
Pharmaceutical Development Services Ltd  
Biocity Nottingham  
Pennyfoot Street  
Nottingham NG1 1GF

This notice has been notified to the appropriate Regulatory Agency.

Wand STA Hand Piece model STA-5050-305.



Wand Hand Piece model WA-2050-305

